



Knockin Stable Cell Lines

Product Manual





Knockin Cell Lines Product Information

1. Description

For knockin stable cell lines, exogenous genes are inserted into specific locations within target cells, which can help researchers study gene functions and establish cell models.

Using our customized service, you will obtain our high-quality customer service experience. The knockin stable cell line is an important research tool for target gene function discovery and research (e.g., signaling pathway research), target validation, antibody research, constructing models, drug screening and assay development. It allows for precise site-directed insertion of exogenous genes and serve as an important tool for studying gene functions and constructing disease models.

2. Knockin Stable Cell Line Generation Service

Service	Technical method	Gene expression evaluation	Application	Delivery
Knockin Stable Cell Line Generation	CRISPR/Cas9	Sequence	Gene function research Protein engineering Disease models creation	Monoclonal cell line

3. Technical advantages


We have optimized and upgraded the CRISPR/Cas9 technology, resulting in the development of a more efficient gene editing system. Our extensive experimental test optimization and cell line construction experiences allow us to enhance the efficiency of homologous recombination during gene knock-in and reduce the probability of random insertion. At present, there have been a large number of successful cases, and more than 1000 gene edits have been successfully achieved in more than 30 kinds of cells. In addition, we provide services such as cellular gene expression regulation and cellular function verification. Finally, professional laboratory reports and quality inspection reports will be delivered.

AcceGen Biotech

277 Fairfield Road,
Fairfield, NJ 07004, USA

T: [+1-862-686-2696](tel:+1-862-686-2696) | E: inquiry@accegen.com
F: [+1-323-978-5598](tel:+1-323-978-5598) | Homepage: www.accegen.com

4. Experimental process



5. Application Operations and Detection Methods

Target gene function discovery and research

- Construct a gene knockin stable cell line
- Culture and expand the knockin stable cell line with the suitable medium and serum
- Group cells according to experimental requirements
- Study the effect of gene knockin on cells through various experimental methods in terms of phenotype or gene function.

Western blotting Cells were lysed using RIPA or other reagents, and the cell supernatant can be taken after centrifugation to determine the concentration of the target protein. Then the difference in protein expression between knockin cells and control cells can be obtained.

RT-PCR The process of extracting cellular RNA can be operated according to the nucleic acid extraction kit, and the target gene product can be obtained after reverse transcription and PCR amplification. Detect and observe the expression results of the target genes by agarose gel electrophoresis and gel imaging system.

Flow Cytometry Inoculate cells (5×10^5 cells/ml) into 6-well plates and culture them for 24 hours. Add antibodies or stimulating factors required for the experiment to the cells and incubate for several hours. Finally, a flow cytometer can be used to detect cell cycle and apoptosis, etc.

Statistical analysis Experimental data can be analyzed using software such as SPSS, GraphPad Prism, Flow Jo, and Excel.

AcceGen Biotech

277 Fairfield Road,
Fairfield, NJ 07004, USA

T: [+1-862-686-2696](tel:+1-862-686-2696) | E: inquiry@accegen.com
F: [+1-323-978-5598](tel:+1-323-978-5598) | Homepage: www.accegen.com



+1-862-686-2696
277 Fairfield Road,
Fairfield, New Jersey 07004, USA

inquiry@accegen.com
www.accegen.com